★ Estimation of alcohol ★

Ethanol (C₂H₅OH) is a primary alcohol that readily undergoes oxidation reactions. The potassium dichromate method is used as a standard analytical method for its determination for several reasons, the most important of which are:

- ❖ High sensitivity in detection.
- ❖ Reliance on a redox reaction.
- ❖ Ability to achieve high-precision quantitative analysis.

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1.Theoretical Background

The reaction is a Redox Reaction where:

- Ethanol → is Oxidized.
- Potassium Dichromate → is Reduced


Substance	Role	Change
Ethanol	Reducing Agent	Changes to Acetic Acid
K ₂ Cr ₂ O ₇	Oxidizing Agent	Changes from Cr ⁶⁺ to Cr ³⁺
Color	Interaction Guide	Orange → Green

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
3.Reaction Chemistry

Role of Acidic Medium

The presence of concentrated H_2SO_4 is essential to:

1. Provide H^+ ions for the reaction.
2. Stabilize the Dichromate ion.
3. Prevent side reactions

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4.Reaction Mechanism .

Why do primary alcohols oxidize so easily:

1. Formation of Chromate Ester between alcohol and Cr_6^{+} .
2. Electron Transfer occurs.
3. Breaking of the C-H bond.
4. Formation of Acetic Acid.

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5. Instrumentation and Reagents

Instrumentation

- Burette (0.01 mL accuracy).
- Volumetric pipette.
- Conical flask.
- Heating source.

Reagents:

- Standard potassium dichromate solution ($K_2Cr_2O_7$)
- Concentrated sulfuric acid (H_2SO_4).
- Ferrous sulfate ($FeSO_4$) as a reducing agent for titration.
- Indicator (e.g., diphenylamine sulfonate) - optional.

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6. Procedure

Step 1: Taking the Sample

Take a specific volume of ethanol to ensure accurate calculations.

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Step 2: Adding Dichromate

Add a calculated excess amount of $K_2Cr_2O_7$. (Reason: To ensure complete oxidation of all ethanol.

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Step 4: Heating

Heat the mixture to speed up the reaction and increase efficiency

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Step 5: Titration

Add Fe^{2+} (Ferrous solution) from the burette. The Fe^{2+} reduces the remaining (excess) dichromate.

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End Point Chemistry

- ❖ Disappearance of orange (Cr^{6+})
- ❖ Appearance of green (Cr^{3+})

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